

Thornton




NATIONAL AERONAUTICS AND SPACE ADMINISTRATION

PRELIMINARY

SMEAT MEDICAL CHECKLIST

(M074, M112-M115, M133, M487, T003)

**EVA/IVA PROCEDURES BRANCH
CREW PROCEDURES DIVISION**



**MANNED SPACECRAFT CENTER
HOUSTON, TEXAS**

MARCH 21, 1972

SMEAT MEDICAL
CHECKLIST

M074, M112-115, M133, M487, T003

MARCH 21, 1972

Prepared By: Tom George
Tom George
SMEAT Experiment Operations
Engineer; Book Manager

Approved By: John McKee
John McKee
Experiment Operations Specialist

Approved By: David C. Schultz
David C. Schultz
Chief, EVA/IVA Procedures Branch

IT IS REQUESTED THAT ANY ORGANIZATION HAVING
COMMENTS, QUESTIONS, OR SUGGESTIONS CONCERN-
ING THIS DOCUMENT CONTACT JOHN MCKEE, EXPE-
RIMENT PROCEDURES SECTION, CG3, BLDG 4,
TELEPHONE 483-3091 OR TOM GEORGE, MISSION
OPERATIONS, MMC BLDG. 1740 NASA BLVD., TELE-
PHONE 333-4150 x141.

DISTRIBUTION OF THIS DOCUMENT IS CONTROLLED
BY J. W. O'NEILL, FLIGHT PLANNING BRANCH,
CREW PROCEDURES DIVISION, TELEPHONE 483-4271.

ACKNOWLEDGEMENT

<u>SECTIONS</u>	<u>NAME</u>	<u>LOCATION</u>
1,2	D.G. SMITH/MMC	333-4150 ext.216 1740 NASA RD. 1
	M.F. Griffin	483-3091
3	N. Vicoli/MMC	333-4150 ext.216 1740 NASA RD.1
	M.F. Griffin	483-3091
4	J. Merchant/MMC	333-4150 ext.228 1740 NASA RD. 1
	L. Ramon	483-4637
5	R. Miller/MMC	333-4150 ext.214 1740 NASA RD. 1
	R. Nute	483-3091

ANY ORGANIZATION HAVING SPECIFIC COMMENTS
REGARDING A PARTICULAR AREA OF RESPONSIBILITY
SHOULD CONTACT THE INDIVIDUAL (S) LISTED
ABOVE.

i
SMEAT
Medical Experiments Checklists
M074, M112-115, M133, M487, T003

LIST OF EFFECTIVE PAGES

Preliminary 3/21/72

Page Number

Change Date

i thru ii	3/21/72
1-1 thru 1-10	3/21/72
2-1 thru TBD	3/21/72
3-1 thru 3-6	3/21/72
4-1 thru TBD	3/21/72
5-1 thru 5-8	3/21/72

ii
CONTENTS

SMEAT MEDICAL CHECKLIST
(M074, M112-M115, M133, M487, T003)

SPECIMEN MASS MEASURING DEVICE
(M074)

EXPERIMENT PREPARATION	PAGE
EXPERIMENT CALIBRATION	1-1
EXPERIMENT PERFORMANCE	1-3
EXPERIMENT STOWAGE	1-7
	1-9

HEMATOLOGY AND IMMUNOLOGY
(M112-M115)

TBS	2-1
-----	-----

SLEEP MONITORING
(M133)

EXPERIMENT PREPARATION	3-1
EXPERIMENT OPERATION	3-3
POST OPERATION ACTIVITIES	3-5

HABITABILITY CREW QUARTERS
(M487)

TBS	4-1
-----	-----

INFLIGHT AEROSOL ANALYSIS
(T003)

T003 - AEROSOL ANALYSIS	5-1
-------------------------	-----

DATE 3/21/72

SPECIMEN MASS MEASUREMENT DEVICE
PREPARATIONM074
PREP

- 613
- 1 cb EXPERIMENTS WMC SMMD - close (up)
(verify)
 - 2 Obtain long 3/16 Allen wrench and handle
from tool kit
 - 3 Open cabinet - verify SMMD control lever
in LOCK position
 - 4 Release tie down sheet - press down
and outward
 - 5 CAM LOCK/CAM UNLOCK - CAM UNLOCK (1/2
turn CCW with Allen wrench)
 - 6 MASS/OFF/TEMP - TEMP
 - 7 RESET - press
 - 8 Temp readout - verify (non-zero reading)
 - 9 MASS/OFF/TEMP - MASS
 - 10 RESET - press
 - 11 Control lever - RELEASE (hold for
10 sec)
 - 12 Control lever - LOCK
 - 13 Digital readout - verify (non-zero reading)
 - 14 MASS/OFF/TEMP - OFF
 - 15 Return Allen wrench and handle to
tool box

DATE 3-21-72

M074
PREP

THIS PAGE INTENTIONALLY LEFT BLANK

DATE 3-21-72

SPECIMEN MASS MEASUREMENT DEVICE -
CALIBRATION

- W742
- 1 Obtain Medical Experiments Log
 - 2 Control lever - LOCK (verify)
 - 3 CAM LOCK/CAM UNLOCK - CAM UNLOCK (verify)
 - 4 MASS/OFF/TEMP - TEMP
 - 5 RESET - press
 - 6 Record reading
 - 7 MASS/OFF/TEMP - MASS
 - 8 RESET - press
 - 9 Control lever - RELEASE (hold for 10 sec)
 - 10 Control lever - LOCK
 - 11 Record reading
 - 12 Repeat steps 8 thru 11 for total of five measurements
 - 13 Compare the readings for repeatability. The span of 4 out of 5 readings should be less than 20 counts
For example: 2.32451
2.32455
2.32442
2.32437
2.32430

If sufficient repeatability of readings was not obtained, repeat the series of measurements

DATE 3-21-72

M074
CALIB

- M074
CALIB
- 14 Release tie down sheet and place 50 gram calibration mass in recess on center of tray
 - 15 Latch tie down sheet and repeat steps 8 thru 13
 - 16 Unlatch tie down sheet and add 100 gram mass on top of 50 gram mass
 - 17 Latch tie down sheet and repeat steps 8 thru 13
 - 18 Unlatch tie down sheet and add 250 gram mass on top of 100 and 50 gram masses
 - 19 Latch tie down sheet and repeat steps 8 thru 13 for total of five measurements
 - 20 Unlatch tie down sheet and remove 250, 100 and 50 gram masses from tray and place in cabinet

CAUTION: DO NOT PLACE MORE
THAN 500 GRAMS ON TRAY

- 21 Place 500 gram mass in center of tray, latch tie down sheet and repeat steps 8 thru 13
- 22 Unlatch tie down sheet and remove 500 gram mass from tray and stow on knurled calib post - then replace 250, 100 and 50 gram masses on post in that order
- 23 Latch tie down sheet and repeat steps 8 thru 13

DATE 3-21-72

- 24 MASS/OFF/TEMP - TEMP
- 25 RESET - press
- 26 Record reading
- 27 MASS/OFF/TEMP - OFF

DATE 3-21-72

THIS PAGE INTENTIONALLY BLANK

DATE 3-21-72

SPECIMEN MASS MEASUREMENT DEVICE -
PERFORMANCE

- 1 CAM LOCK/CAM UNLOCK - CAM UNLOCK
(Verify)
- 2 Control lever - LOCK (Verify)
- 3 Release tie down sheet
- 4 Place specimen against tray, relatch
tie down sheet
- 5 MASS/OFF/TEMP - TEMP
- 6 RESET - press
- 7 Record reading
- 8 MASS/OFF/TEMP - MASS
- 9 RESET - press
- 10 Control lever - RELEASE (hold until
display stops counting)
- 11 Control lever - LOCK
- 12 Record reading
- 13 Repeat steps 8 thru 11 for total of
five readings

DATE 3-21-72M074
PERFORMANCE

- 14 Compare readings for repeatability
The span of 4 out of 5 readings
should be less than 100 counts
For example:
6.32451
6.32539 - (highest reading)
6.32482
6.32492
6.32429 - (lowest reading)

Note this series exhibits acceptable
repeatability in 4 out of 5 readings

- 15 MASS/OFF/TEMP - TEMP
16 RESET - press
17 Record reading
18 MASS/OFF/TEMP - TEMP
19 Control lever - LOCK (Verify)
20 Obtain wet wipe, unlatch tie down
sheet, clean tray and sheet, relatch
sheet

SPECIMEN MASS MEASUREMENT DEVICE -
STOWAGE

- E624 1 Obtain long 3/16 Allen wrench and
 handle from tool box
- 2 Open cabinet - verify control lever
 in LOCK position
- 3 CAM LOCK/CAM UNLOCK - CAM LOCK (1/2
 turn CW with Allen wrench)
- E624 4 Return long 3/16 Allen wrench and
 handle to tool box

DATE 3-21-72M074
STOWAGE

THIS PAGE INTENTIONALLY BLANK

DATE 3-21-72

M074
STOWAGE

M112 - M115

TBS

DATE 3/21/72

THIS PAGE INTENTIONALLY LEFT BLANK

DATE 3/21/72

EXPERIMENT EQUIPMENT PREPARATION

- D418 Obtain high power cable
- 551 Connect high power cable to HI POWER
ACCESSORY OUTLET, 5 or 6
- Thread high power cable thru floor grid
and penetration in sleep station 1
light baffle
- F578 Release three fasteners holding M133
panel launch stowage bracket top
frame
- F578 Remove M133 panel assembly from launch
stowage bracket
- Secure stowage bracket top frame to base
- Move M133 panel assembly to sleep station 1
- 904 cb EXP PWR - open (down) (verify)
- ELECTRODE SELECT - OFF (verify)
- TAPE RECORDER - 1 (verify)
- SLEEP STA 1 Attach M133 panel assembly to dovetail
fixture face toward bulkhead
- S913 Obtain M133 power cable and SIA cable
- Connect power cable to J5 and SIA cable
to J3 on panel assembly
- E624 Using Phillips screwdriver, engage
recorder pinch roller on tape recorder
1 (outer recorder)

DATE 3/21/72

M133 PREP

SLEEP Remove M133 panel assembly from dove-
STA 1 tail fixture, reverse and install
(face outward)

Connect high power cable to M133 power
cable

901 Connect SIA cable to CHAN B connector
on SIA

904 cb EXP PWR - CLOSE (up)

M133 PREP

DATE _____

EXPERIMENT OPERATION

S913 Obtain cap, preamp-accel assembly and chin strap.

SLEEP Using scissors, remove tips from elec-
STA 1 trodes. Dispose of tips and cap bag in trash container

S913 Connect preamp-accel assembly to cap

SLEEP Don cap, secure using chin strap
STA 1

Enter sleep restraint

Connect preamp-accel cable to CAP UMB receptacle on panel

904 SUBJECT GAIN - predetermined position

cb EXP PWR - close (up) (verify)

ELECTRODE SELECT - TEST

ELECTRODE STATUS (6 indicators) - lighted

Note: If either EEG channel can not be lighted by rocking the appropriate electrodes obtain new cap

DATE 3/21/72

M133 OPERATION

904 ELECTRODE SELECT - LEFT EEG

TAPE RECORDER -1 (verify)

The subject may begin sleep period

Note: If necessary to exit sleep station before sleep period is complete, turn ELECTRODE SELECT - OFF and disconnect preamp-accel cable from panel. Upon return, reconnect cable and check electrode status before resuming recording and sleep.

The following steps are performed upon awaking.

ELECTRODE SELECT - TEST

ELECTRODE STATUS (6 indicators) -
lighted

Note: If any ELECTRODE STATUS indicators are not lighted, note in experiment logbook.

ELECTRODE SELECT - OFF

SLEEP Remove cap
STA 1

Remove preamp-accel assembly from cap
and disconnect preamp-accl cable
from the panel

S913 Place preamp-accl assembly and chin
strap in stowage cabinet

S190 Dispose of cap

SLEEP Using wipes, remove any electrolyte
STA 1 remaining in hair

DATE 3/21/72

POST OPERATION ACTIVITIES

904 cb EXP PWR - open (down)

SLEEP Remove panel assembly from mount
STA 1

Reverse panel assembly and replace on
dovetail fixture (face toward bulkhead)

Remove outer recorder from panel assembly

E624 Obtain #2 Phillips screwdriver

SLEEP Using Phillips screwdriver, disengage
STA 1 the pinch roller

Using Phillips screwdriver, loosen six
screws securing recorder cover and
remove cover

Using Phillips screwdriver, press reel
hubs and remove both reels

S913 Obtain tape return canister

SLEEP Remove new reels from tape return
STA 1 canister

Place old reels in tape return canister

Note: Use care not to deform or twist
tape. Refer to recorder cover
diagram for supply and takeup
reel locations and tape path.

DATE 3-21-72

MT33 POST OPER

SLEEP Using Phillips screwdriver to press the
STA 1 reel hub, install supply and takeup
reels in recorder

Thread tape thru recorder tape path

Place cover on recorder and secure six
retaining screws

Replace tape recorder on panel
assembly

Remove panel assembly from mount

Reverse panel assembly and replace on
dovetail mount (face outward)

4-1

M487

M487

TBS

DATE 3/21/72

M487

THIS PAGE INTENTIONALLY LEFT BLANK

DATE 3/21/72

T003 AEROSOL ANALYSIS

DAILY MEASUREMENT

- 1 Remove AEROSOL ANALYZER from locker 613
Replace data card if card has less than 7
blanks
- 2 Translate to EXPERIMENT COMP near wardroom
wall
- 3 Set filter - impactor dial to position #1
- 4 Record GMT and filter position
Hold INLET port near ceiling and
perpendicular to longitudinal axis

INITIATE CYCLE button - push

Observe PILOT OVERFLOW light - off after 70
seconds

Record CHANNEL COUNT - #1
Record CHANNEL COUNT - #2
Record CHANNEL COUNT - #3

- 5 Stow AEROSOL ANALYZER IN LOCKER 613

DATE 3-21-72

T003 AEROSOL
ANALYSIS DAILY
MEASUREMENT

T003 AEROSOL
ANALYSIS DAILY
MEASUREMENT

THIS PAGE INTENTIONALLY LEFT BLANK

DATE 3-21-72

10 DAY MEASUREMENTS

- 1 Remove AEROSOL ANALYZER from locker 613
Replace data card if card has less than 7
blanks
- 2 Translate to EXPERIMENT COMP near wardroom
wall
- 3 Set filter - impactor dial to position #1
- 4 Record GMT and filter position
Hold INLET port near ceiling and
perpendicular to longitudinal axis

INITIATE CYCLE button - push

Observe PILOT OVERFLOW light - off after 70
seconds

Record CHANNEL COUNT - #1

Record CHANNEL COUNT - #2

Record CHANNEL COUNT - #3

- 5 Translate to EXPERIMENT COMP near LBNPD
Set filter - impactor dial to position #2
Repeat step 4
- 6 Translate to EXPERIMENT COMP near MANLOCK
Set filter - impactor dial to position #3
Repeat step 4
- 7 Translate to SLEEP COMP
Set filter - impactor dial to position #4
Repeat step 4
- 8 Translate to WASTE MANAGEMENT AREA
Set filter - impactor dial to position #5
Repeat step 4

DATE 3-21-72

10 DAY
MEASUREMENT

- 9 Translate to WARDPOOM near table
Set filter - impactor dial to position #6
Repeat step 4
- 10 Translate to UPPFP DECK AREA
Set filter - impactor dial to position #7
Repeat step 4
- 11 Stow AEPOSOL ANALYZER in locker 613

Note: Filter position #8 is for ASTPONAUTS
discretion (20 measurements max)

10 DAY
MEASUREMENT

DATE 3-21-72

MEAL MEASUREMENT

- 1 Remove AEROSOL ANALYZER from locker 613
Replace data card if card has less than 2
blanks
- 2 Translate to WARDROOM near table
- 3 Set filter - impactor dial to position #6
- 4 Record GMT and filter position
Hold INLET port near ceiling and
perpendicular to longitudinal axis

INITIATE CYCLE button - push

Observe PILOT OVERFLOW light - off after 70
seconds

Record CHANNEL COUNT - #1

Record CHANNEL COUNT - #2

Record CHANNEL COUNT - #3

- 5 Stow AEROSOL ANALYZER nearby universal
camera bracket
After meal repeat step 4
- 6 Stow AEROSOL ANALYZER in locker 613

DATE 3-21-72

MEAL
MEASUREMENT

THIS PAGE INTENTIONALLY LEFT BLANK

MEAL
MEASUREMENT

DATE 3-21-72

CLOTHING MEASUREMENT

- 1 Remove AEROSOL ANALYZER from locker 613
Replace data card if card has less than 2
blanks
- 2 Translate to area where clothes are to be
changed
- 3 Set filter - impactor dial to position #4
- 4 Record GMT and filter position
Hold INLET port near ceiling and
perpendicular to longitudinal axis

INITIATE CYCLE button - push

Observe PILOT OVERFLOW light - off after 70
seconds

Record CHANNEL COUNT - #1
Record CHANNEL COUNT - #2
Record CHANNEL COUNT - #3

- 5 Change clothes
- 6 Repeat step 4

Stow AEROSOL ANALYZER in locker 613

DATE 3-21-72

THIS PAGE INTENTIONALLY LEFT BLANK

NASA MSC

DATE 3-21-72