JUNE 18, 1974

5.2.3 M074, Specimen Mass Measurement

Functional Objectives - The objectives of experiment M074 were to demonstrate mass measurement accuracy without gravity, to validate theoretical behavior of the device to measure fluids as well as solids, and to support M071/M073 experiments requiring mass determination.

Experiment Description

The Specimen Mass Measurement Device (SMMD) used a linear spring pendulum. Measurement of mass was based on the inertial properties of mass. The device measured and recorded the time associated with the period of the plate - fulcra spring - supported pendulum that has a fixed displacement.

The restrained mass to be measured was accelerated uniformly by a repeatable restoring force set into the plate-fulcra springs. Three periods of a pendulum were timed. The spring period was measured electro-optically. The measurement was displayed electronically and was converted to a direct mass readout.

- Concept The SMMD incorporated a linear spring mass pendulum to measure mass based on the inertial properties of mass. This device measured and recorded the time associated with the period of a plate-fulcra spring supported pendulum with a fixed displacement. The masses to be measured were uniformly accelerated by a repeatable restoring force set into the plate-fulcra springs, and three periods of a pendulum were timed. Known masses were measured and a nomograph developed which could be used to determine the mass of other objects. This device had a minimum response to gravitational force which was eliminated in orbital flight by recalibration with the same known masses as used preflight.
- 5.2.3.3 <u>Key Personnel</u> Principal Coordinating Scientist:

Whittle, M.W., M.D. Squadron Leader, Royal Air Force

Principal Investigator:

Thornton, W.E., M.D., Johnson Space Center, Houston, Texas

5.2.3.4 Experiment Performance Data - The M074 experiment consisted of calibrating the SMMD with specified weight periodically during each Skylab mission as follows:

	<u>SL-2</u>	<u>SL-3</u>	<u>SL-4</u>
Inflight Calibration	5	5	6
Special Tests		2	

5.2.3.4.1 Preflight Data

A A CO

- A. Formal KSC Checkout This checkout verified SMMD operational status performed final pre launch calibration.
- B. Preflight Calibration Data During preflight the SMMD was calibrated a number of times with a standard set of masses. These masses were later used in fligt calibration.

5.2.3.4.2 Inflight Data

- A. Calibration Data The SMMD was calibrated inflight with the same standard set of masses used in preflight calibrations. These measurements were used to provide a corrective factor for the gravitational force effects on the preflight measurements.
- B. Food Residue Data Food not completely consumed by a crewman during the mission was measured with SMMD and recorded.
- C. Fecal Data Each fecal deposit with wipes was measured with the SMMD and recorded during the mission.
- D. Special Tests Tests were performed to measure the SMMD repeatability.
- E. Crew Comments SMMD operation crew comments were recorded on the dump tapes as part of the M487 experiment.

5.2.3.4.3 Postflight Data

A. Crew Comments - SMMD operation crew comments and performance are part of the Crew Technical Debriefings for each mission.

5.2.3.4.4 Postflight Reports

- A. Principal Investigator's Special M074 Studies In preparation.
- B. Principal Investigator's MO74 Final Report